

Fermod F430 Latch Set - No Lock & 20mm Packer. For Door thickness 67mm - 120mm.

PART No: F430CRH-120P

FEATURES

Features:

- European design & styling
- Antimicrobial Resistant
- Independently Tested and shown to reduce Corona Viruses & Bacteria
- Constructed from Industrial Grade Composite for strength plus negative & positive temperature conditions.
- Invisible & Tamperproof fixing
- Interchangeable for left & right-hand opening doors
- Safe and easily activated Internal release mechanism
- Soft & quiet operation



Combine this latch set with Fermod F1573 Hinges & a Kason K109308 Door Closer for the perfect self-closing door operation.

CRH & All manufacturers highly recommend that all internal safety release mechanisms fitted to doors are checked at least 3 times weekly to ensure they are functioning correctly.

Please view the video tab below for a demonstration



Colour: Grey & Black

Material: Composite

ADDITIONAL INFORMATION

All Fermod latch sets now include an antimicrobial resistant treatment on the external & internal lever that limits the growth of harmful microorganisms*.

Tests conducted by independent laboratories according to ISO 21702 (for coronavirus**) and ISO 22196 (for bacteria) how that fasteners incorporating the feature allow:

- An 87% reduction in the number of coronaviruses after 1 hour and 96% after 8 hours
- A 99% reduction in the number of bacteria after 1 hour
- As the feature is an integral part of the material, the product is protected over its lifetime, without loss of efficiency.
- The treatment is compatible with food contact and in testing shown to

reduce the growth of bacteria responsible for food poisoning on the surface. These fasteners are perfectly in line with the HACCP process.

- More effective than a simple anti-bacterial treatment, because it acts on coronavirus, fungi and mould.
- Totally odorless and invisible, the addition of antimicrobial chemicals does not affect the quality, integrity or function of the fasteners.
- This treatment has been used for many years and is perfectly safe.
- The use of antimicrobial resistant fasteners complements good practices and barriers to limit the spread of viruses and bacteria and cannot replace them.
- The tests were carried out on the HCoV-229E.
- Reduction of the treated sample against the control of 75% after 1 hour and 92% after 8 hours.

Set consists of Latch Body, Internal Release Handle, Rods, Fasteners & Strike Plate.



CRH Australia Pty Ltd 1-3 Nelson Street Moorabbin VIC 3189 Australia t+61 3 9532 f+61 3 9532 2696

info@crh.com.au

www.crh.com.au